

Product datasheet for TP521118

Asphd1 (NM_001039645) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse aspartate beta-hydroxylase domain containing 1 (Asphd1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR221118 representing NM_001039645 Red =Cloning site Green =Tags(s)
	MSANTFGNAGFSVLLPGARLEGRCGPTNARVRCHLGLKIPPGCELVGGEPQCWAEGHCLLVDDSLHTV AHNGSPEDGPRVVFIVDLWHPNVAGAERQALDFVFAPDP
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	11.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001034734
Locus ID:	233879
UniProt ID:	Q2TA57
RefSeq Size:	518
Cytogenetics:	7 F3
RefSeq ORF:	327



[View online »](#)

Synonyms: A830007L07Rik; Gm168