

Product datasheet for **TP521304**

Pnmal2 (NM_001099636) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse PNMA-like 2 (Pnmal2), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR221304 representing NM_001099636 Red =Cloning site Green =Tags(s)
	 MTMNULLQDWCKELEVEVHRALLITGIPERLQQADIEATLSPNLQPLGNYRLRTVRAVTKEKSQAALVEFG ADFDHSAIPIHIRGENGIWKILSKDRGQDARVLRQMRRLLLDERPMDDFREIAIPLVLEAQAQAKGLGKE AKKAGSSEGAARSGRRGRRGGKRRSRSHKLTAQKGQKRSRGGRRARSESEDSSDDSLGIVIEEIIYAEDLS VDEDQRALYATLQAAAKELTKKWAFRGDQEEGDGPREFLALVTVTDKAKKEIEKDLPGTESICLNKDE KSGVPDLVALLAVRDTPAVEVETEVDDEDEDEDEDEDEDEDEDEDEDESDSELENGEHGKEVMDSP FVAIVAYTDPADPSAREEMLKIASVIETLWGDGKDKKDVLPQVLSVMAKDTSGPRVKVEEAGRQVDAMV LRKAEEDGNLLECISTLAEAENCPKGKKSGLGLLRGWTAEGHQGGLLELVALLAAQDMVEAVREEEASRW GSGGSGGRKCDHSQGGLEVLAFSLASQENLESNEESEDSDTESEGDSEDTDSESEPGDSVSKKPRAKR ARTGSKGLPPAGATAVSTASRARKARRGGRGRGRAVTPEKKTGSGVSTEDHAGNNNNKKKKSAGAGARAR AGEAKGQPAAVPKSTRGKKARRGPRRAPKCR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	71.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_001093106](#)
Locus ID: 434128
UniProt ID: [G3X9N3](#)
RefSeq Size: 4147
Cytogenetics: 7 A2
RefSeq ORF: 1983
Synonyms: EG434128