

Product datasheet for TP521323

Morn1 (NM_001081100) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse MORN repeat containing 1 (Morn1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR221323 representing NM_001081100 Red =Cloning site Green =Tags(s) MEAAGDPYQASLQPRRHSPRLPPRNGYGVVYPNSFFRYEGEWKGGKKHGHGKLLFKDGSYYEGEFVDG E ITGEGCQHWAWSGNTYSGQFVLGEPQGHGVMKYKAGGRYEGELSQGLREGQGFLEDQDGQVYQGSFH DNK RHGRGQMIFKNGDKYEGDWVRDQRQGHGVLLCCADGSTYKGQWHSDFVSGLSLVHCSGVTYCGMFIN GHP AEQAKKIVILGPEVLEVQGSPTLSVQLQQDDGEVAKSESGRVLKISAGVRYVQLPEYSEVSFFKMDDED YMETPIQTFFGFQCISYPLSVPMWSWGPEPGSTVESARADLLLSKRDSEPLDLSGAFHGGQDGLSILPARR HKPCCQQVDHGCDEFVDILLGAPPPGMQPYLFLPGTLEKTLGRP KGGPSSPEVLPTAQEPLKGTNRPDGT TVEPPTAAYLGEYVLMVCDVTSPFLGHRLPTTFKHLRLAKGDTHWPHVPEDHQKL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	52.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq: NP_001074569
Locus ID: 76866
UniProt ID: A2RTS7
RefSeq Size: 1661
Cytogenetics: 4 E2
RefSeq ORF: 1431
Synonyms: 2900057D20Rik; 4930417P05Rik; AU042805