

## Product datasheet for TP521347

### Uba7 (NM\_023738) Mouse Recombinant Protein

#### Product data:

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Mouse ubiquitin-like modifier activating enzyme 7 (Uba7), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>MR221347 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MDEELYSRQLYVLGLPAMQRIQEAKVLLCGLQQLGAEVAKNLVLTGVGSLTLHDPHPTCWADLAAQCFLS EESLGRNRAEASQAQLAQLNEAVQISVHSGDITEDLLQGFQVWVLTDSKLEDQLKVGPLCHKHGVRFLMA ETRGLVGRFLCDFGEDFTVLDPTEVEPMATAIQDISQGFPGIVTLRGDTKRHSFHDGDLVIFSDIEGMVE LNSCSPQSVRVQKDGSLGIDTTTFSRYLRGGVVTEVKRPKTVRHKPLDIALQLPHVVAQNTQEVQRAHC LHQAFHVLHKFQQLHGRLPKPWPDDAETVELAQDLEPLKGTTEEESLDEALLRTIALSSAGTLPMAAI MGGVAAQEVKKAISRKFMPLDQWLYFDALECLPEDETLLPSPEDCQPRNCRYDGGQIAVFGTDLQEKLSQD HYLLVGAGAIGCEMLKVFLVGLGVRANGGVTADMDYIERSNLSRQFLFRPKDVRPKAEVAAAAAHL NPDLRATPYTCPLDPTTEDIYDDSFVRVNGVVAALDSFQARHYVAARCTHYLKPILLEAGTQGTWGSASV FVPYVTEAYRGPASDAASEDAPYPVCTLRHIPSSMEHSVQWAQDQFEGFLRSTETINCYQQTCTLSAT DRTETLALLQQVMGVLTRPQTWQDCVWALGHWQLCFHDKVLEGGTQFSSGSNKCPHPLQFDPNHDMHF LYVLAANLYARMHGLPGSQSPALRELLTRLES SRPQNLFSAEHGQEQLKELQETLDDWRKGPPLKP VLFVKDDDSNFHVDVVAATDLRCQNYGILPVNHARIKQIVGRIIPAIATSTAVVAGLLGLELYKVVSGL RSHGTFRHSYLHLAENHFIRSAPSAPAVQSFRDLKWTWCWDRLLKVPVAVQPRTLKSLLAHLQEEHGLKVEM LLHHQALLYSSGWSSEKQAQHLCLRVTELQHVTVGWKPKPGLKVLVLFELSCEGEEEEEMAFPLHYEL</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
<b>Tag:</b>	C-MYC/DDK
<b>Predicted MW:</b>	108.7 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C after receiving vials.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_076227</a>
<b>Locus ID:</b>	74153
<b>UniProt ID:</b>	<a href="#">Q9DBK7</a>
<b>RefSeq Size:</b>	3162
<b>Cytogenetics:</b>	9 59.07 cM
<b>RefSeq ORF:</b>	2934
<b>Synonyms:</b>	1300004C08Rik; Ube1l