

## Product datasheet for TP521490

### Sec14I3 (NM\_001029937) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse SEC14-like lipid binding 3 (Sec14I3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR221490 representing NM_001029937 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MSGRVGDLSPKQAETLAKFRENVQDVLPALPNPDDYFLLRWLRARNFDLQKSEAMLRKYMFRKTMDID H ILDWQPPEVIQKYMPPGGLCGYDRDGCPVWYDIIGPLDPKGLLFSVTKQDLLKTKMRDCERILHECDLQTE RLGRKIETIVMIFDCEGLGLKHFVKPLVEVYQEFFGLLEENYPETLKFMLIVKATKLPVGYNLMKPFLS EDTRRKIVVLGNSWKEGLLKLISPEELPAHFGGTLTDPDGNPKCLTKINYGGEIPKSMYVRDQVKTQYE HSVQISRGSSSHQVEYEILFPGCVLRWQFSSDGADIGFGVFLKTKMGERQKAGEMVEVLTSQRYNAHMPVE DGS LTCSEAGVYVLRFDNTYSFVHAKKVSFTVEVLLPDEGMQKYDEELTPV
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	46.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_001025108</a></u>
Locus ID:	380683



[View online »](#)

UniProt ID: [Q5SQ27](#)  
RefSeq Size: 3928  
Cytogenetics: 11 A1  
RefSeq ORF: 1203  
Synonyms: 1110069O07Rik