

Product datasheet for TP521548

Tex13a (NM_026469) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse testis expressed 13A (Tex13a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR221548 representing NM_026469 Red=Cloning site Green=Tags(s) MALNPEDLGSGFRHSKVSFMFINEQMAKHDKGPDFYLENLSLSWEEVEDKFKVILEDSEVSNEAREACAW G TLALGVRFAHRQGCLQGHGVQWLQDLSSLHKMSALSLSPLDKQLIHQQEMAAQKEVALQLQMAQAKLE EVQ RERDLLRLKILQAE LRALPNAVRPAVAIPPAVVRGGIRTQWSSTKENLAEWMAAAATGRPNERANMSDTA LSGTITSPEEVLEDPNGSFMKLLGVMKCKNYPLKRQLDLRPKKASMCSFSQSLNLRSTVSSEPFIVQLP ASFTYSYESPFPAIPTTSQLPTTERQPHVYPYSMASDISHLSDMRIHRGDHQLPKDKRFSAFRRPGDWD CPWCKAVNFSRRENCFHCGKGIWLQNP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	42.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_080745</u>


[View online »](#)

Locus ID: 67944

UniProt ID: [A2AFS9](#)

RefSeq Size: 1317

Cytogenetics: X F1

RefSeq ORF: 1131

Synonyms: 1700025D03Rik; RP23-149B17.1