

## Product datasheet for TP522328

### H1f1 (NM\_030609) Mouse Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse histone cluster 1, H1a (Hist1h1a), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR222328 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MSETAPVAQAASTATEKPAAAKKTKKPAKAAAPRKKPAGPSVSELIVQAVSSSKERSGVSLAALKKSLAA AGYDVEKNNSRIKLGKSLVNGTLVQTKGTGAAGSFKLNKKAESKAITTKVSVKAKASGAAKPKKTAG AAAKKTVKTPKKPKPAVSKKTSKSPKPKVVKAKKVAKSPAKAKAVKPKASKAKVTKPKTPAKPKKAAP KKK
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-MYC/DDK
Predicted MW:	21.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_085112</a>
Locus ID:	80838
UniProt ID:	<a href="#">P43275</a>
RefSeq Size:	745



[View online »](#)

**Cytogenetics:** 13 A3.1

**RefSeq ORF:** 642

**Synonyms:** H1; H1-1; H1.1; H1a; H1v; H1var3; Hist1; Hist1h1a

**Summary:** Histones are basic nuclear proteins responsible for nucleosome structure of the chromosomal fiber in eukaryotes. Two molecules of each of the four core histones (H2A, H2B, H3, and H4) form an octamer, around which approximately 146 bp of DNA is wrapped in repeating units, called nucleosomes. The linker histone, H1, interacts with linker DNA between nucleosomes and functions in the compaction of chromatin into higher order structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H1 family. Transcripts from this gene lack polyA tails but instead contain a palindromic termination element. [provided by RefSeq, Aug 2015]