

Product datasheet for TP522434

Triml1 (NM_177742) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse tripartite motif family-like 1 (Triml1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR222434 representing NM_177742 Red =Cloning site Green =Tags(s) MSNHEKMSTTDLMENLREELTCFICLDYFSSPVTECGHSFCLMCLLKSWEEHNTPLSCPECWRTLGA FQANERLGRLANIGRQLRSQVLQSEDEQSICGRMPGPSWVFSDDDEQSVINVSPPSQGTNKACFSSEAE EQHKEKLQDIINILRKKKKEVQAILNHEKERVMLCKEETKTCKQVVVSEYMKMHQFLKEEQQLQLL EREEKANMKKLRENEIQLTQQIRRLGKMIGRIESTCQNLTLSEFEEVKGAMDRIESLLFQSPETTIT ELSLCHITGMREMLRKFDITLDPATANAYLLLEDLKSVMRYGGTRQQLPDNPFRDQSATVLGAQIF TCGRHYWEVEVGKTEWEVGICKDSVNRKGNLPKPPGDLFSLIGLKIGDDYSLWVSSPLKGQHVR EPHVKGVFLDYD SGHIAFYFNATDESLIYFPPTPFHEALRPIFSPCLPNEGNTNDPLIICH TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	53.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_808410
Locus ID:	244448


[View online »](#)

UniProt ID: [Q8BVP1](#)

RefSeq Size: 1850

Cytogenetics: 8 A4

RefSeq ORF: 1410

Synonyms: 4933403D14; BC050188

Summary: Probable E3 ubiquitin-protein ligase which plays an important role in blastocyst development. Involved in progression of blastocyst stage and subsequent embryo development.[UniProtKB/Swiss-Prot Function]