

Product datasheet for TP523161

Tcp11l1 (NM_177190) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse t-complex 11 like 1 (Tcp11l1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR223161 protein sequence Red =Cloning site Green =Tags(s) MSENLDKSHVDEAGEAEEAASEQGLEGALECSDETLQKKVKSDSPSSQRVGRPHSSPARLVTVEELLETA KGVTNMLAHEIVVTGDFRINAVELAEGSLEKRVKEIVHKAFWDCLSVQLSEEPPTYDHAIKLVGEIKET LLSFLLPGHTRLRNQITEVLDLELIKQEAENGALDISKLAEFIIGMMGILCAPARDEEVKKLKGIKEIVP LFRAIFSVLDLMKVDMANFAISSIRPHLMQQSVEYERRKFQEVLERQPNSLDFATQWLEEATNDLLSQKY KHALPAGGGAAGSGDAPLLTPVSVQNFAYLKLLKWDHFQRPFPETVLMDDQSRFQELQLQLEQLAVLGAV L LVTFTAAPGISGHADFAEKLKMMVKTLTDMHLPFSFHLEDALASIGEKVCLEVSSCLSLCGSSPFSVAK ETVLKGQIQALASPEDPIRRIMESRIFTLETYLASGHQKPLPTVPGGLGPIQKELEEAVKVFRLVNYN KMFVFCPPYYDAILSKLLRS TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	56.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq:	<u>NP_796164</u>
Locus ID:	320554
UniProt ID:	<u>Q8BTG3</u>
RefSeq Size:	3512
Cytogenetics:	2 E2
RefSeq ORF:	1527
Synonyms:	C130096D04Rik