

Product datasheet for **TP524102**

Fnbp4 (NM_018828) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse formin binding protein 4 (Fnbp4), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR224102 representing NM_018828 Red =Cloning site Green =Tags(s)
	MRKCTAPVKSGSPREAVPRFMNIDYAPGRCAPAPLICRRSLLALRLMMGKKSRAVPGRRPILQLSPPGPR SSTPGRDPDPDPDEADSTAAATSQSAPAAATAAAATSPAVPASAAPEDSPSEDEQEWWVEVPNVVQNPP TPVMTTRPTAVKATGGLCLLGAYADSDDDESVDSEKTAQSKEANGNQATDIDSTLANFLAEIDAITAPQP AAPVVASAPPPTPPRPEPKEAATPALSPTASNGSDTAQTPGWHYDTQCSLAGVEIEMGDWQEWWDENTG C YYYWNTQTNEVTWELPQYLATQVQGLQHYQPSSVTGTEAAFVNTDMYTKERTTAASSSKSGPVITKREV KKEVNEGIQALSNEEERKGVAAALLAPLLPEGVKEEEERWRRKVICKADPVSETKETSTASEETGPSI KPPEVMMDGTEDPSQEELCSVVQSGESEEEEEEEQDTLELELALERKKAELRALEEGDGSVSGSSPRSD ISQPASQDGVRRIMSKRGKWKMFVRATSPSTSRSSSKTGRDSPENGETAIGAEDSEKIDEISDKETEVE ESSEKIKVQLAPKVEEQDLKFQIGELANTLSKFEFLGINRQSSNFHMLLLQTETRIADWREGALNGN YLKRKLQDAAEQLKQYEINATPKGWSCHWDRDHRRYFYVNEQSGESQWEFPDGEESQTKVRDESL P KLTVKDKTCTDPNSTESSENPTGSLCKESFSGQVSSSLMPLTPFWTLLQSNVPVLQPPLPLEMPPPPPP PESPPPPPPPPPPPLEDGEIQEVEMEDEGSEEPAPGTEEDTPLKPSTQTTAVTSQSLVDSTASSPPS NKAVKRKAPEMSTSVVQRSATIGSSPVLYSQSAIAAGHQAVGMAHQAVGMAHQAVSASHAAAAGVGHQ AR GMSLQSNYLGAAAPALMSYAECVPIGVTTPSLQPAQARGTMAAPAVVEPPPPPPPPPTPTPPPPPPAP KVPPPEKTRKGKKDKAKKSKTKMPSLVKKWQSIQRELDEEDNSSSEEDRESTAQKRIEWWKQQQLVSGL AERNANFEALPEDWRARLKRRKMAPST TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	116.3 kDa



[View online »](#)

Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_061298
Locus ID:	55935
UniProt ID:	Q6ZQ03
RefSeq Size:	4691
Cytogenetics:	2 E1
RefSeq ORF:	3231
Synonyms:	FBP30; Fnbp30; mKIAA1014