

Product datasheet for **TP524164**

Jrk (NM_008415) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse jerky (Jrk), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR224164 representing NM_008415 Red =Cloning site Green =Tags(s)

MASKQAAAKGKGEKRKRVLTLKEKIDICTRLERGESRKALMQEYNVGMSTLYDIKAHKAQLLRFASDD
SRQALEQRRTLHTPKLEHLDRVLYEWFLVKRAEGIPVSGPMLIEKAKDFYKQMRLEPCVFSGGWLWRFK
ARHGIKKLDASSEKQAADHQAAEQFCGFFRSLAAEHGLSPEQVYSADETGLVWRCLPNSAPDDGTVPHFK
QGKDRLTVLMCANATGSHRIKPLAIGKGGGPRAFRGIQHLPIAYKAQGNAWVDKEIFSDWFHHIFVPSVR
EHFRTIGLPEDSKAILLLDHSRAHSQESELVSENVFTIFLPSSVTSLLQPTEQGIRRAFMRLEFINPPVAF
QGFPTRHINDAIVNVARAWNAVPSQVFQRAWRKLWPTVTFTEGSSEEEAECCAIPKHTFAHILGLVK
EGPTCSGSRQLQDSRVEERVVAGRDIDEAPAIVAPSQATRCKKSEKDTGETEEAAWEQAATSFEALVRFA
ERQPCFSVQEMGQLQALHTVFRRQQQLRQPRVALRAVIKLEALQEHPGVCVATHTPLPCSSTAGDN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	62.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_032441



[View online »](#)

Locus ID: 16469

UniProt ID: [Q60976](#), [Q543Z9](#)

RefSeq Size: 3378

Cytogenetics: 15 34.26 cM

RefSeq ORF: 1671

Summary: May bind DNA.[UniProtKB/Swiss-Prot Function]