

Product datasheet for **TP524253**

Hspbap1 (NM_175111) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse Hspb associated protein 1 (Hspbap1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR224253 representing NM_175111 Red =Cloning site Green =Tags(s)
	MAAGCEGIAPPTLGERTVGEEGEPVKPFTPEKAKEIIMSLQQPAIFCNMVFWDWPSRHWTAKHLSKVLEGK QIRFRMGLRGTVGPQYETEC SYVDATLEEF LTWNC DQSSISGPFKDYEH SKFWAYADYKYFVTLFEDKT DVFQEVWSDFGFPGRNGQESTLWIGSFGAHTPCHLDSYGCNLV FQVQGRKRWHLFPPEDTPFLYPTRI P YEESVFSKINVNPDLKCFPQFQKARRHMTLS PGQVLFVPRHWWHYVESLDPVTVSINSWIELEDHL ARVEEAITRMLVCTLKTAEDPHHPRTWLNPTVEEETSHEVNSCYLNSAVCAFFDHCEKAKAVELQVLSAN GAEP RVQEEHMEVEQAQGPSLTSRAEKQEA VSLFGPDLPVTPASEERGGALEGDSEESVSSNGGHFAEL PCARRQQTSGASALTEPLAPRGPMASRVFVSTDDLLDCLVNPQVTRMVAQLLIQ GKSCDTF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	54.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_780320</u>



[View online »](#)

Locus ID:	66667
UniProt ID:	<u>Q8BK58</u>
RefSeq Size:	4606
Cytogenetics:	16 B3
RefSeq ORF:	1449
Synonyms:	3830421G21 Rik
Summary:	May play a role in cellular stress response.[UniProtKB/Swiss-Prot Function]