

Product datasheet for TP524272

Klhdc1 (NM_178253) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse kelch domain containing 1 (Klhdc1), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR224272 representing NM_178253 Red=Cloning site Green=Tags(s) MADPQAFCAEERSGHC AVVDGHFLYVWGGYVSIEDNEVYLPNDEMWTYDIDSGLWKMHLMEGELPPS MS GSCGACIHGRLYVFGGYDDKGYSNRLYFVNLRTDGTYTWEKITKFDGQPPTPRDKLSCWVYKDRLIYFG GYGYRRHSELQECFDVHDASWEEQIFWGW HNDVHVFDTKTRTWSQPEIKGGVPPQPRAAHSCAVLGNK GY VFGGRVLQTRMNDLHYLNLDTWWSGRISVNGESPKHRSWHTLTAITDDKLFLFGGLNADNIPLSDGWI H NITTNCWKQLRHLPTYPRLWHTACLGKENEIMVFGGSKDNLLFLDTGHCNDLLIFQTQPYSLLRSCLDL IGKNAILKLSQISLLPPKLLQVLKKITFWTAANYRKEQRIKKEETENNQPRVSSC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	46.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.


[View online »](#)

RefSeq:	<u>NP_839984</u>
Locus ID:	271005
UniProt ID:	<u>Q80YG3</u>
RefSeq Size:	1949
Cytogenetics:	12 C2
RefSeq ORF:	1218
Synonyms:	MGC141300; MGC141301