

Product datasheet for TP524357

Cbx3 (NM_007624) Mouse Recombinant Protein

Product data:

| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Purified recombinant protein of Mouse chromobox 3 (Cbx3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug |
| Species: | Mouse |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >MR224357 representing NM_007624 Red=Cloning site Green=Tags(s) |

MASNKTTLQKMGKKQNGKSKKVEEAPEEFVVEKVLDRRVNGKVEYFLKWKGFADADNTWEPEENLDCP
ELIEAFLNSQKAGKEKDGTKRKSLSDESDDSKSKKRDAADKPRGFARGLDPERIIGATDSSGELMFLM
KWKDSDEADLVLAKEANMKCPQIVIAFYEERLTHWSCPEDIAQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

| | |
|----------------|--|
| Tag: | C-MYC/DDK |
| Predicted MW: | 20.8 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C after receiving vials. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | <u>NP_031650</u> |
| Locus ID: | 12417 |
| UniProt ID: | <u>Q9DCC5</u> , <u>Q32P00</u> |
| RefSeq Size: | 1851 |
| Cytogenetics: | 6 24.89 cM |



[View online »](#)

RefSeq ORF: 549

Synonyms: HP1g; M32