

Product datasheet for TP524359

Trim47 (NM_172570) Mouse Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse tripartite motif-containing 47 (Trim47), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR224359 representing NM_172570 Red =Cloning site Green =Tags(s) MDGSGPFSCPICLEPLREPVTLP CGHNFCLACLGALWPHRSAGGTGGSGGPARCPLCQEPFPDGLQLRK N HTLSELLQLRQGSVPGPMSAPASGSTRGATPEPSAPSAPPPAPEPSAPCAPEQWPAGEEPPVRCDACPEGA ALPAALSCLSLCLASFCSAHLAPHERSPALRGHRLVPPLRRLEESLCPRHLRPLERYCRVERVCLCEACAT QDHRGHELVPLEQERALQEVEQSKVLSAAEDRMDELGAGIAQSRRTVALIKSAAVAERERSQMFAEATA TLQSFQNEVMGFIEEGEATMLGRSQGDLRREQEQRSLSKARHNLGQVPEADSVSFLQELLALRLALEEG CGPGPGPPRELSFTKSSQVVKAVRDTLISACASQWEQLRGLGSNEDGLQKLGSEDVESQDPDSTSLLESE APRDYFLKFAYIVDLSDTADKFLQLFGTKGVKRVLCPINYPESPTRFTHCEQVLGEGALDRGTYWVEVE IIEGWVSVGVMAEGFSPQEPYDRGRLGRNAHSCCLQWNGRGFSVWFCGLEAPLPHAFSPTVGCLEYAD H ALAFYAVRDGKLSLLRRLKASRRSGALASPTDPFQSRLDSHFSGLFNHLKPAFFLESVDAHLQIGPL KKSCITVLKRRSGP SGP TRRRLE EQKLISEEDLA AND I LDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	70.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.


[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_766158</u>
Locus ID:	217333
UniProt ID:	<u>Q8C0E3</u>
RefSeq Size:	2192
Cytogenetics:	11 E2
RefSeq ORF:	1932
Synonyms:	2210023F24Rik
Summary:	E3 ubiquitin-protein ligase that mediates the ubiquitination and proteasomal degradation of CYLD.[UniProtKB/Swiss-Prot Function]

Product images:



Purified recombinant protein Trim47 was analyzed by SDS-PAGE gel and Coomassie Blue Staining.