

Product datasheet for TP524602

Armc3 (NM_001081083) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse armadillo repeat containing 3 (Armc3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR224602 representing NM_001081083 Red =Cloning site Green =Tags(s) MGKKIKKEVEPPPKDVFDPITIESKKAATVVLMLKSPEEDILAKACEAIYKFALKGEENKATLLELGAVE PLTKLLTHEDKTVRRNAMMIFGILASNSDVKKLLRELEVMSVIAQLSPEEEVVIHEFASLCLANMSVEY TGKVQIFEHGGLEPLIRLLSSSDPDVKNSIECIYNLVQDFQCRTTLQELNAIPPILELLRSEYPIIQLL ALKTLGVITCDKEARTMLKENQGLDHLTKILETKELNDLHVEALSIVANCLEDMDTMVLMQQTGSLKKVL SFAESSTIPDIQNAAKAIKAAAYDPENRKVFHEQEVEKCLVTLLGSDSDGTKIAASQAISALCENLSCK EFFNTQGIPQIVQLLRSDNEEVREAAALANLTTSSPANANAAAEADAIDPLINILSSKRDGAIANAAT VLTNMATQEPLRAIQNHMHALGPLHSTNTLVQSTAALTVAATACDVEARTQLRNCGGLVPLVGLLH SKNDEVRRHASWAVMVCAGDEPMAVELCRLGALNILEEINRSLSRKNKFSEAAAYNKLLNNNLSLKYSQTG YLSSSNIISDGFYDYGRINPGTKLLSLKELCLQELNDQRAILLVNNKSDTSPPPSMEDKSSDVGYGRSIS SSSSLRRGSKEKANAIFGSPTEEKSEPASVRNTILSRAFTKEKGVSKKKSRLQLICSSYLLWKGKGKKEE EKVKEEEEILALPKFTEGSPEKEWNPPDPPEFCVYVLEVTKSILPIVNLKEQIEVLAKYVADKMGGKIPK EKLADFSWELHISELKFQLKSNVPIGYIKKGIFYHRALLFKALADKIGVGCSLVRGEYSRGWNEVKLVN EARKGMIGNLPPPEEYIVDLMFHPGNLLKLSKEADLYRFL TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-MYC/DDK
Predicted MW:	97.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


[View online »](#)

Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001074552</u>
Locus ID:	70882
UniProt ID:	<u>A2AU72</u>
RefSeq Size:	3004
Cytogenetics:	2 A3
RefSeq ORF:	2643
Synonyms:	4921513G22Rik