

## Product datasheet for TP526093

### Nfatc3 (NM\_010901) Mouse Recombinant Protein

#### Product data:

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Purified recombinant protein of Mouse nuclear factor of activated T cells, cytoplasmic, calcineurin dependent 3 (Nfatc3), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
<b>Species:</b>	Mouse
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>MR226093 representing NM_010901 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MTTANCGAHDELDFKLVFGEDGAPAPVSQVSRPADLEPDDCASIFNVDPPPSTLNSSLGLPHHGLLQS HSSVLSPSFQLQGYKNEYGTGDISISKYSPLGGPKPFECPSIQITSISPNCHQGTDAHEDDLHINDPERE YLERPSRDHLYLPLEPSYRESSLSPSPASSISSRSWFSDASSCESLSHIYDDVDSELNEAAARFTLGSP TSPGGSPGGCPGEESWHQQYSGHSLSPRQSPCHSPRSSITDENWLSRPASGPPSRPTSPCGKRRHSSA EVCYAGSLSPHHSPVSPGHSPRGSVTEDTWLTAPVHTGSGLSAPFPFYCVETDIPLKTRKTSQAA ILPGKLEICSDDQGNLSPSRETSVDDGLGSQYPLKKDSSGDQFLSVSPFTWSKPKPGHTPIFRTSSLPP LDWPLPTHFGQCELKIEVQPKTHHRAHYETEGSRGAVKASTGGHPVVKLLGYSEKPINLQMFIGHTADDRY LRPHAFYQVHRITGKTVATASQEIISTKVLKLEIPLPENNMASIDCAGILKLRNSDIELRKGETDIGR KNTRVRLVFRVHIPQPSGKVLSLQIASIPVECSQRSQELPHIEKYSINSCSVNGGHEMIVTGSNFLPES KIIFLEKQDGRPHWEAEGKIIEKQGAHIVLEVPPYHNPAVTSAVQVHFYLCNGKRKKSQSQRFTYTP VLMKQEQRQEDTDLPSVSLPVPVHSAQAQRPSSETGHPHDRAMSAAGSLLCQVQPAYTSMVASTHLPQLQC RDEGAGKEQHIATSSVMHQPFQVTPTSPIGSSYQSIQTSMYNGPTCLPVNPASSQEFDPVLFQQDAALSS LVNLGCQPLSPIPFHSSNSDATGHLLAHSVSVQTPPHLQSMGYHCSNAGQALSSPVADQITGQPSSHL QPITYCPSHPGSATAASPAASHPLASSPISGPPSQLPMPYQSPSSGTASSPSPTRMHSGQHSTQAQS TGQGGLSVPSLVCHSLCDPASFPFGGATVSIKPEPEDQEPNFATIGLQDITLDDVNEIIGRDMSSQISVS QATEVMRDTPLPGPASPDLMTSHSAQ</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
<b>Tag:</b>	C-MYC/DDK
<b>Predicted MW:</b>	116 kDa
<b>Concentration:</b>	>0.05 µg/µL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining



[View online »](#)

<b>Buffer:</b>	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C after receiving vials.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_035031</a>
<b>Locus ID:</b>	18021
<b>UniProt ID:</b>	<a href="#">Q3UZ64</a>
<b>RefSeq Size:</b>	6032
<b>Cytogenetics:</b>	8 53.08 cM
<b>RefSeq ORF:</b>	3228
<b>Synonyms:</b>	C80703; D8Erttd281e; NFAT4; NFATx