

Product datasheet for TP526421

Ddx20 (NM_017397) Mouse Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 20 (Ddx20), with C-terminal MYC/DDK tag, expressed in HEK293T cells, 20ug
Species:	Mouse
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>MR226421 representing NM_017397 Red=Cloning site Green=Tags(s)

MAAAAFEVPAALTTSESTMAAERAAAPVQAVEPTASPWTQRTAHDIGGPRTRTGDVWLAEPADFESLLL
SRPVLEGLRAAGFERPSPVQLKAIPLGRCGLDLIVQAKSGTGKTCVFSTIALDSLILENYSTQILILAPT
REIAVQIHSVITAIGIKMEGLECHVFIGGTPLSQDKTRLKCKHIAVGSPIGRIQLIELDYLNPGSIRLFI
LDEADKLLLEEGSFQEINWIYSSLPASKQMLAVSATYPEVLANALTRYMRDPTFVRLNPSDPSLIGLKQY
YQVNSYPLAHKIFEKTQHLQELFSKVPFNQALVFSNLHSRAQHLADILSSKGFPTTECISGNMNQNRQL
DAMAKLKQFHCRLVLISTDLTSRGIDAEKVNLDVPLDWETYMHRIGRAGRFGTLGLTVTYCCRGEEE
NMMMIAQKCNINLLPLDPPIPPGLMEECLNWDVEVKAAMHTYSSPTVATQSPKKQVQKLERAFQSQRTP
GNQTPSPRNTSASALSARPKHSPKLPVKSHSECGVLEKAAPPQESGCPAQLEEQVKNVQTSVEDSSSN
SQHQAKDSSPGSLPKIPCLSSFKVHQPSTLTFAELVDDYEHYIKEGLEKPEIIRHYTGPEAQTGNPQNG
FVRNRVSEDRAQMLVSSSQSGDSESDSCSSRTSSQSKGNKSYLEGSSDTQLKDTECTPVGGPLSLEQV
QNGNDTPTQVEYQEAPETQVKARHKEGANQRSKQSRNPARRSSYRVQSEPQEEWYDCHRETTASFSDT
YQDYEEYWRAYYRAWQEYAAAASHSYWNAQRHPSWMAAYHMNTVYLQEMMRGNQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-MYC/DDK
Predicted MW:	92.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C after receiving vials.



[View online >](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_059093
Locus ID:	53975
UniProt ID:	Q9JY4
RefSeq Size:	2710
Cytogenetics:	3 F2.2
RefSeq ORF:	2475
Synonyms:	dp103; GEMIN3
Summary:	<p>The SMN complex plays a catalyst role in the assembly of small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome. Thereby, plays an important role in the splicing of cellular pre-mRNAs. Most spliceosomal snRNPs contain a common set of Sm proteins SNRPB, SNRPD1, SNRPD2, SNRPD3, SNRPE, SNRPF and SNRPG that assemble in a heptameric protein ring on the Sm site of the small nuclear RNA to form the core snRNP. In the cytosol, the Sm proteins SNRPD1, SNRPD2, SNRPE, SNRPF and SNRPG are trapped in an inactive 6S pICln-Sm complex by the chaperone CLNS1A that controls the assembly of the core snRNP. Dissociation by the SMN complex of CLNS1A from the trapped Sm proteins and their transfer to an SMN-Sm complex triggers the assembly of core snRNPs and their transport to the nucleus. May also play a role in the metabolism of small nucleolar ribonucleoprotein (snoRNPs) (By similarity).[UniProtKB/Swiss-Prot Function]</p>