

## Product datasheet for **TP700140**

### Eph receptor A3 (EPHA3) (NM\_182644) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of human EPH receptor A3 (EPHA3), transcript variant 2, with C-terminal DDK/His tag, expressed in human cells, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from TrueORF clone, RC600040, encoding the region (Glu21- Val539) of human EPHA3
Tag:	C-DDK/His
Predicted MW:	62 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH 7.4, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_872585</a>
Locus ID:	2042
UniProt ID:	<a href="#">P29320</a>
RefSeq Size:	2684
Cytogenetics:	3p11.1
RefSeq ORF:	1617
Synonyms:	EK4; ETK; ETK1; HEK; HEK4; TYRO4



[View online »](#)

**Summary:**

This gene belongs to the ephrin receptor subfamily of the protein-tyrosine kinase family. EPH and EPH-related receptors have been implicated in mediating developmental events, particularly in the nervous system. Receptors in the EPH subfamily typically have a single kinase domain and an extracellular region containing a Cys-rich domain and 2 fibronectin type III repeats. The ephrin receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands. This gene encodes a protein that binds ephrin-A ligands. Two alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:**

Druggable Genome, Protein Kinase, Secreted Protein, Transmembrane

**Protein Pathways:**

Axon guidance

**Product images:**