

## Product datasheet for **TP700143**

### Eph receptor A6 (EPHA6) (NM\_001080448) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of human EPH receptor A6 (EPHA6), transcript variant 1, with C-terminal DDK/His tag, expressed in human cells, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from TrueORF clone, RC600043, encoding the region (Ala23 – Ser549) of human EPHA6
Tag:	C-DDK/His
Predicted MW:	62 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH 7.4, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001073917</a>
Locus ID:	285220
UniProt ID:	<a href="#">A0A0B4J1T8</a> , <a href="#">B3KS12</a>
RefSeq Size:	3677
Cytogenetics:	3q11.2
RefSeq ORF:	1647
Synonyms:	EHK-2; EHK2; EK12; EPA6; HEK12; PRO57066



[View online »](#)

**Summary:**

Receptor tyrosine kinase which binds promiscuously GPI-anchored ephrin-A family ligands residing on adjacent cells, leading to contact-dependent bidirectional signaling into neighboring cells. The signaling pathway downstream of the receptor is referred to as forward signaling while the signaling pathway downstream of the ephrin ligand is referred to as reverse signaling (By similarity).[UniProtKB/Swiss-Prot Function]

**Protein Families:**

Druggable Genome, Protein Kinase, Transmembrane

**Protein Pathways:**

Axon guidance

**Product images:**