

## Product datasheet for **TP700164**

### CCK4 (PTK7) (NM\_152880) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of human protein tyrosine kinase 7 (PTK7), transcript variant 2, with C-terminal DDK/His tag, expressed in human cells, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from TrueORF clone, RC600064, encoding the region (Ala31-Thr704) of human PTK7
Tag:	C-DDK/His
Predicted MW:	73 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH 7.4, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_690619</a>
Locus ID:	5754
UniProt ID:	<a href="#">Q13308</a>
RefSeq Size:	4153
Cytogenetics:	6p21.1
RefSeq ORF:	1992
Synonyms:	CCK-4; CCK4



[View online »](#)

**Summary:**

This gene encodes a member of the receptor protein tyrosine kinase family of proteins that transduce extracellular signals across the cell membrane. The encoded protein lacks detectable catalytic tyrosine kinase activity, is involved in the Wnt signaling pathway and plays a role in multiple cellular processes including polarity and adhesion. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]

**Protein Families:**

Druggable Genome, Protein Kinase, Transmembrane

**Product images:**