

Product datasheet for **TP700197**

CD47 (NM_001777) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens CD47 molecule (CD47), residues 19-141 aa, with C-terminal DDK/His tag, expressed in HEK293 cells.
Species:	Human
Expression Host:	HEK293
Tag:	C-DDK/His
Predicted MW:	35.21
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH7.4, 10% glycerol
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001768
Locus ID:	961
RefSeq Size:	5346
Cytogenetics:	3q13.12
RefSeq ORF:	969
Synonyms:	IAP; MER6; OA3
Summary:	This gene encodes a membrane protein, which is involved in the increase in intracellular calcium concentration that occurs upon cell adhesion to extracellular matrix. The encoded protein is also a receptor for the C-terminal cell binding domain of thrombospondin, and it may play a role in membrane transport and signal transduction. This gene has broad tissue distribution, and is reduced in expression on Rh erythrocytes. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2010]
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	ECM-receptor interaction



[View online »](#)

Product images:

