

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TP700221

Lymphocyte Activation Gene 3 (LAG3) (NM_002286) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of human lymphocyte-activation gene 3 (LAG3), with C-terminal DDK/His tag, expressed in human cells, 20 μg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	A DNA sequence from TrueORF clone, RC220269, encoding the region (Val29-Leu450) of human LAG3
Tag:	C-Fc
Predicted MW:	70 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS, pH 7.4, 10% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP 002277</u>
Locus ID:	3902
UniProt ID:	<u>P18627</u>
RefSeq Size:	1991
Cytogenetics:	12p13.31
RefSeq ORF:	1575
Synonyms:	CD223




This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	Lymphocyte Activation Gene 3 (LAG3) (NM_002286) Human Recombinant Protein – TP700221
Summary:	Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like

mmary:Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like
domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and
chromosomal localization all indicate a close relationship of LAG3 to CD4. [provided by
RefSeq, Jul 2008]

Protein Families: Transmembrane

Product images:



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US