


## Product datasheet for TP701124

### SARS-CoV-2 Spike-RBD Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant SARS coronavirus 2 Spike-RBD Protein (Arg319~Phe541), with C-terminal Mouse FC tag,secretory expressed in HEK293 cells, 50ug
Species:	SARS-CoV-2
Expression Host:	HEK293
Tag:	C-Fc Mouse
Predicted MW:	51.3 KDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">QHD43416</a></u>

#### Product images:



[View online »](#)