


Product datasheet for TP701125

SARS-CoV-2 Spike-RBD Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant SARS coronavirus 2 Spike-RBD protein(Lys310~Leu560) with C-terminal His tag,secretory expressed in HEK293 cells, 50ug
Species:	SARS-CoV-2
Expression Host:	HEK293
Tag:	C-HIS
Predicted MW:	30.1 KDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	PBS
Storage:	Store at -80°C after receiving vials.
Stability:	Stable for at least 3 months from receipt of products under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	QHD43416

Product images:



Purified recombinant protein S Protein (RBD) was analyzed by SDS-PAGE gel and Coomassie Blue Staining.



[View online »](#)