

## Product datasheet for TP701169

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## SARS-CoV-2 Spike-RBD Protein (A372S) Mutant MN908947) COVID-19 Recombinant Protein

**Product data:** 

**Product Type: Mutant Proteins** 

Description: Purified recombinant SARS coronavirus 2 Spike-RBD protein(Lys310~Leu560),A372S, with C-

terminal His tag, secretory expressed in HEK293 cells, 50ug

Species: SARS-CoV-2

**Expression Host:** HEK293T

**Expression cDNA Clone** 

KGIYQTSNFR VQPTESIVRF PNITNLCPFG EVFNATRFAS VYAWNRKRIS NCVADYSVLY or AA Sequence: NSSSFSTFKC YGVSPTKLND LCFTNVYADS FVIRGDEVRQ IAPGQTGKIA DYNYKLPDDF

TGCVIAWNSN NLDSKVGGNY NYLYRLFRKS NLKPFERDIS TEIYQAGSTP CNGVEGFNCY FPLQSYGFQP TNGVGYQPYR VVVLSFELLH APATVCGPKK STNLVKNKCV NFNFNGLTGT

**GVLTESNKKF L** 

C-HIS Tag:

Predicted MW: 30.1 KDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 95% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** PBS, 10 mM PB, pH 7.2, 150 mM NaCl

Note: For culture applications, please filter before use. Note that you may experience some loss of

protein during the filtration process.

Storage: Store at -80°C after receiving vials.

Stable for 12 months from the date of receipt of the product under proper storage and Stability:

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: QHD43416