


Product datasheet for **TP710279**

HHV-6Au96 Human herpesvirus Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of HHV-6Au94, with C-terminal DDK tag, expressed in sf9, 20ug
Species:	Human herpesvirus
Expression Host:	Sf9
Tag:	C-DDK
Predicted MW:	55.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	50 mM Tris-HCl, 100 mM glycine, pH 8.0, 150 mM NaCl, 30% glycerol
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

Product images:



[View online >](#)