

Product datasheet for TP720033XL

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

Il1b (NM 008361) Mouse Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of mouse Interleukin-1B/IL-1B

Species: Mouse Expression Host: E. coli

Expression cDNA Clone

Val118-Ser269

or AA Sequence:

Tag: Tag Free
Predicted MW: 17.4 kDa
Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Bioactivity: ED50 is less than 0.01 ng/ml as determined by the dose-dependent stimulation of mouse

D10S cells. Specific Activity of 1.0 x 108 IU/mg.

Endotoxin: < 0.1 EU per μg protein as determined by LAL test

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

RefSeg: NP 032387

 Locus ID:
 16176

 UniProt ID:
 P10749

 RefSeq Size:
 1328

Cytogenetics: 2 62.97 cM

RefSeq ORF: 807

Synonyms: IL-; Il-1b; IL-1beta





Summary:

The protein encoded by this gene is a member of the interleukin 1 cytokine family. This cytokine is produced by activated macrophages as a proprotein, which is proteolytically processed to its active form by caspase 1. The encoded protein plays a role in thymocyte proliferation and is involved in the inflammatory response. [provided by RefSeq, Aug 2015]

Product images:

