

Product datasheet for TP720035

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Fgf2 (NM_008006) Mouse Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of mouse fibroblast growth factor basic

Species: Mouse Expression Host: E. coli

Expression cDNA Clone

Met1-Ser154

or AA Sequence:

Met 1-Set 154

Tag: Tag Free
Predicted MW: 17.1 kDa
Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per μg protein as determined by LAL test

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 μ g/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 6 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 032032

Locus ID: 14173

UniProt ID: <u>P15655</u>, <u>Q541T2</u>

RefSeq Size: 695

Cytogenetics: 3 18.41 cM

RefSeq ORF: 462

Synonyms: bFGF; Fgf-2; Fgfb





Summary:

Acts as a ligand for FGFR1, FGFR2, FGFR3 and FGFR4. Also acts as an integrin ligand which is required for FGF2 signaling. Binds to integrin ITGAV:ITGB3. Plays an important role in the regulation of cell survival, cell division, cell differentiation and cell migration. Functions as a potent mitogen in vitro. Can induce angiogenesis.[UniProtKB/Swiss-Prot Function]

Product images:

