

Product datasheet for **TP720092XL**

CA7 (NM_001014435) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human carbonic anhydrase VII (CA7), transcript variant 2
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Ala264
Tag:	C-His
Predicted MW:	30.7 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0 .
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_001014435
Locus ID:	766
UniProt ID:	P43166
Cytogenetics:	16q22.1
Synonyms:	CA-VII; CAVII
Summary:	Carbonic anhydrases are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization. The cytosolic protein encoded by this gene is predominantly expressed in the brain and contributes to bicarbonate driven GABAergic neuron excitation. Alternative splicing in the coding region results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2018]



[View online »](#)

Protein Families: Druggable Genome

Protein Pathways: Nitrogen metabolism

Product images:

