

Product datasheet for **TP720097XL**

Lon Escherichia coli Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of E. coli Lon Protease
Species:	Escherichia coli
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Asn2-Lys784
Tag:	C-His
Predicted MW:	88.5 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 μ m filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.



[View online »](#)