

Product datasheet for **TP720128XL**

LACTB2 (NM_016027) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human lactamase, beta 2 (LACTB2)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Leu288
Tag:	N-GST
Predicted MW:	59.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Supplied as a 0.2 um filtered solution of PBS,1mMDTT PH7.5.
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at < -20°C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_057111
Locus ID:	51110
UniProt ID:	Q53H82 , A0A024R811
Cytogenetics:	8q13.3
Synonyms:	CGI-83
Summary:	Endoribonuclease; cleaves preferentially 3' to purine-pyrimidine dinucleotide motifs in single-stranded RNA. The cleavage product contains a free 3' -OH group. Has no activity with double-stranded RNA or DNA. Required for normal mitochondrial function and cell viability. [UniProtKB/Swiss-Prot Function]



[View online »](#)

Product images:

