

## Product datasheet for **TP720132**

### Apoptosis repressor with CARD (NOL3) (NM\_001185057) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nucleolar protein 3 (apoptosis repressor with CARD domain) (NOL3), transcript variant 3.
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Ser208
Tag:	Tag Free
Predicted MW:	22.6 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001171986</a>
Locus ID:	8996
UniProt ID:	<a href="#">O60936</a> , <a href="#">O60936-1</a>
Cytogenetics:	16q22.1
Synonyms:	ARC; FCM; MYOCL1; MYP; NOP; NOP30
Summary:	This gene encodes an anti-apoptotic protein that has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2010]



[View online »](#)

**Product images:**

