

## Product datasheet for **TP720133XL**

### **BNP (NPPB) (NM\_002521) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human natriuretic peptide precursor B (NPPB)
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	His27-His134
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	14.2 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
<b>Endotoxin:</b>	< 0.1 EU per µg protein as determined by LAL test
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
<b>RefSeq:</b>	<a href="#">NP_002512</a>
<b>Locus ID:</b>	4879
<b>UniProt ID:</b>	<a href="#">P16860</a>
<b>Cytogenetics:</b>	1p36.22
<b>Synonyms:</b>	BNP; Iso-ANP



[View online »](#)

**Summary:**

This gene is a member of the natriuretic peptide family and encodes a secreted protein which functions as a cardiac hormone. The protein undergoes two cleavage events, one within the cell and a second after secretion into the blood. The protein's biological actions include natriuresis, diuresis, vasorelaxation, inhibition of renin and aldosterone secretion, and a key role in cardiovascular homeostasis. A high concentration of this protein in the bloodstream is indicative of heart failure. The presence of myocardial injury is a significant predictor of mortality in hospitalized coronavirus disease 2019 (COVID-19) patients, and there is evidence of increased levels of natriuretic peptide B in hospitalized non-survivor COVID-19 patients. The protein also acts as an antimicrobial peptide with antibacterial and antifungal activity. Mutations in this gene have been associated with postmenopausal osteoporosis. [provided by RefSeq, Aug 2020]

**Protein Families:**

Druggable Genome, Secreted Protein, Stem cell - Pluripotency

**Product images:**