

## Product datasheet for **TP720137XL**

### DNA Polymerase beta (POLB) (NM\_002690) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human polymerase (DNA directed), beta (POLB)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ser2-Glu335
Tag:	C-His
Predicted MW:	39.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 $\mu$ m filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per $\mu$ g protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_002681</a>
Locus ID:	5423
UniProt ID:	<a href="#">P06746</a>
Cytogenetics:	8p11.21
Summary:	The protein encoded by this gene is a DNA polymerase involved in base excision and repair, also called gap-filling DNA synthesis. The encoded protein, acting as a monomer, is normally found in the cytoplasm, but it translocates to the nucleus upon DNA damage. Several transcript variants of this gene exist, but the full-length nature of only one has been described to date. [provided by RefSeq, Sep 2011]
Protein Families:	Druggable Genome
Protein Pathways:	Base excision repair



[View online »](#)

**Product images:**

