

Product datasheet for TP720138

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PPM1G (NM 177983) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human protein phosphatase 1G (formerly 2C), magnesium-

dependent, gamma isoform (PPM1G), transcript variant 1

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met317-Asp546

Tag: C-His

Predicted MW: 26.1 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per µg protein as determined by LAL test

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 817092

Locus ID: 5496

UniProt ID: <u>O15355</u>, <u>Q6IAU5</u>

Cytogenetics: 2p23.3

Synonyms: PP2CG; PP2CGAMMA; PPP2CG

Summary: The protein encoded by this gene is a member of the PP2C family of Ser/Thr protein

phosphatases. PP2C family members are known to be negative regulators of cell stress response pathways. This phosphatase is found to be responsible for the dephosphorylation of Pre-mRNA splicing factors, which is important for the formation of functional spliceosome. Studies of a similar gene in mice suggested a role of this phosphatase in regulating cell cycle

progression. [provided by RefSeq, Apr 2010]





Protein Families: Druggable Genome, Phosphatase

Product images:

