

Product datasheet for TP720140

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PPP1R14A (NM_033256) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human protein phosphatase 1, regulatory (inhibitor) subunit 14A

(PPP1R14A)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met1-Pro147

Tag: C-His

Predicted MW: 17.8 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per µg protein as determined by LAL test

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 150281

 Locus ID:
 94274

 UniProt ID:
 Q96A00

 Cytogenetics:
 19q13.2

Synonyms: CPI-17; CPI17; PPP1INL

Summary: The protein encoded by this gene belongs to the protein phosphatase 1 (PP1) inhibitor family.

This protein is an inhibitor of smooth muscle myosin phosphatase, and has higher inhibitory activity when phosphorylated. Inhibition of myosin phosphatase leads to increased myosin phosphorylation and enhanced smooth muscle contraction. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene. [provided by RefSeq, Sep

20111





Protein Families: Druggable Genome

Protein Pathways: Vascular smooth muscle contraction

Product images:

