

Product datasheet for TP720165L

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PGP9.5 (UCHL1) (NM_004181) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human ubiquitin carboxyl-terminal esterase L1 (ubiquitin

thiolesterase) (UCHL1)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met1-Ala223

Tag: C-His

Predicted MW: 25.9 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per µg protein as determined by LAL test

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 004172

Locus ID: 7345

UniProt ID: <u>P09936</u>, <u>V9HW74</u>

Cytogenetics: 4p13

Synonyms: HEL-117; HEL-S-53; NDGOA; PARK5; PGP 9.5; PGP9.5; PGP95; SPG79; Uch-L1

Summary: The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol

protease that hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in the neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease.[provided by RefSeq, Sep

2009]

Protein Families: Druggable Genome, Protease





Protein Pathways: Parkinson's disease

Product images:

