

## Product datasheet for **TP720166L**

### UCHL3 (NM\_006002) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ubiquitin carboxyl-terminal esterase L3 (ubiquitin thiolesterase) (UCHL3)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Met1-Ala230
Tag:	C-His
Predicted MW:	27.2 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_005993</a>
Locus ID:	7347
UniProt ID:	<a href="#">P15374</a> , <a href="#">A0A140VJZ4</a>
Cytogenetics:	13q22.2
Synonyms:	UCH-L3
Summary:	The protein encoded by this gene is a member of the deubiquitinating enzyme family. Members of this family are proteases that catalyze the removal of ubiquitin from polypeptides and are divided into five classes, depending on the mechanism of catalysis. This protein may hydrolyze the ubiquitinyl-N-epsilon amide bond of ubiquitinated proteins to regenerate ubiquitin for another catalytic cycle. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Aug 2012]



[View online »](#)

Protein Families: Druggable Genome, Protease

**Product images:**

