

## Product datasheet for **TP720182L**

### Annexin V (ANXA5) (NM\_001154) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human annexin A5 (ANXA5)
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ala2-Asp320
Tag:	Tag Free
Predicted MW:	35.9 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_001145</a>
Locus ID:	308
UniProt ID:	<a href="#">P08758</a> , <a href="#">V9HWEO</a>
Cytogenetics:	4q27
Synonyms:	ANX5; ENX2; HEL-S-7; PP4; RPRGL3



[View online »](#)

**Summary:**

The Annexin 5 gene spans 29 kb containing 13 exons, and encodes a single transcript of approximately 1.6 kb and a protein product with a molecular weight of about 35 kDa. The protein encoded by this gene belongs to the annexin family of calcium-dependent phospholipid binding proteins some of which have been implicated in membrane-related events along exocytotic and endocytotic pathways. Annexin 5 is a phospholipase A2 and protein kinase C inhibitory protein with calcium channel activity and a potential role in cellular signal transduction, inflammation, growth and differentiation. Annexin 5 has also been described as placental anticoagulant protein I, vascular anticoagulant-alpha, endonexin II, lipocortin V, placental protein 4 and anchorin CII. Polymorphisms in this gene have been implicated in various obstetric complications. [provided by RefSeq, Dec 2019]

**Product images:**