

## Product datasheet for **TP720201L**

### GFAP (NM\_002055) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens glial fibrillary acidic protein (GFAP), transcript variant 1
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Leu292-Met432
Tag:	N-His
Predicted MW:	18.7 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_002046</a>
Locus ID:	2670
UniProt ID:	<a href="#">P14136</a> , <a href="#">A7REI1</a>
Cytogenetics:	17q21.31
Synonyms:	ALXDRD
Summary:	This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008]
Protein Families:	ES Cell Differentiation/IPS



[View online »](#)

**Product images:**

