

## Product datasheet for **TP720202XL**

### **A2LD1 (GGACT) (NM\_033110) Human Recombinant Protein**

#### **Product data:**

<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human AIG2-like domain 1 (A2LD1)
<b>Species:</b>	Human
<b>Expression Host:</b>	E. coli
<b>Expression cDNA Clone or AA Sequence:</b>	Met1-Arg153
<b>Tag:</b>	N-His
<b>Predicted MW:</b>	19.5 kDa
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>95% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
<b>Endotoxin:</b>	< 0.1 EU per µg protein as determined by LAL test
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
<b>RefSeq:</b>	<a href="#">NP_149101</a>
<b>Locus ID:</b>	87769
<b>UniProt ID:</b>	<a href="#">Q9BVM4</a>
<b>Cytogenetics:</b>	13q32.3
<b>Synonyms:</b>	A2LD1
<b>Summary:</b>	The protein encoded by this gene aids in the proteolytic degradation of crosslinked fibrin by breaking down isodipeptide L-gamma-glutamyl-L-epsilon-lysine, a byproduct of fibrin degradation. The reaction catalyzed by the encoded gamma-glutamylaminocyclotransferase produces 5-oxo-L-proline and a free alkylamine. Two transcript variants encoding the same protein have been found for this gene.[provided by RefSeq, Aug 2010]



[View online »](#)

Product images:

