

Product datasheet for **TP720217XL**

PAFAH1B2 (NM_001184746) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human platelet-activating factor acetylhydrolase 1b, catalytic subunit 2 (30kDa) (PAFAH1B2), transcript variant 2.
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ser2-Ala229
Tag:	C-His
Predicted MW:	26.6 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_001171675
Locus ID:	5049
UniProt ID:	P68402
Cytogenetics:	11q23.3
Synonyms:	HEL-S-303
Summary:	Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene. [provided by RefSeq, Jan 2014]
Protein Pathways:	Ether lipid metabolism, Metabolic pathways



[View online »](#)

Product images:

