

## Product datasheet for TP720223L

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## PMM2 (NM\_000303) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human phosphomannomutase 2 (PMM2)

Species: Human
Expression Host: E. coli

**Expression cDNA Clone** 

Met1-Ser246

or AA Sequence:

Tag: C-His

Predicted MW: 29.1 kDa

Concentration: lot specific

**Purity:** >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

**Endotoxin:** < 0.1 EU per μg protein as determined by LAL test

**Storage:** Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeg: NP 000294

**Locus ID:** 5373

UniProt ID: <u>015305</u>, <u>A0A0S2Z4J6</u>, <u>Q59F02</u>

Cytogenetics: 16p13.2

Synonyms: CDG1; CDG1a; CDGS; PMI; PMI1; PMM 2

Summary: The protein encoded by this gene catalyzes the isomerization of mannose 6-phosphate to

mannose 1-phosphate, which is a precursor to GDP-mannose necessary for the synthesis of dolichol-P-oligosaccharides. Mutations in this gene have been shown to cause defects in glycoprotein biosynthesis, which manifests as carbohydrate-deficient glycoprotein syndrome

type I. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome





**Protein Pathways:** 

Amino sugar and nucleotide sugar metabolism, Fructose and mannose metabolism, Metabolic pathways

## **Product images:**

