

Product datasheet for TP720228

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CUTA (NM 001014433) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human cutA divalent cation tolerance homolog (E. coli) (CUTA),

transcript variant 1

Species: Human
Expression Host: E. coli

Expression cDNA Clone

or AA Sequence:

Met1-Pro156

Tag: C-His

Predicted MW: 17.9 kDa **Concentration:** lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per µg protein as determined by LAL test

Reconstitution Method: Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the

lyophilized protein in ddH2O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeq: NP 001014433

 Locus ID:
 51596

 UniProt ID:
 060888

 Cytogenetics:
 6p21.32

Synonyms: ACHAP; C6orf82

Summary: May form part of a complex of membrane proteins attached to acetylcholinesterase (AChE).

[UniProtKB/Swiss-Prot Function]

Protein Families: Transmembrane





Product images:

