

Product datasheet for **TP720242L**

NADK (NM_001198994) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NAD kinase (NADK), transcript variant 3.
Species:	Human
Expression Host:	E. coli
Expression cDNA Clone or AA Sequence:	Ser64-Gly446
Tag:	N-His
Predicted MW:	44.4 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Storage:	Store at -80°C.
Stability:	Stable for at least 3 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_001185923
Locus ID:	65220
UniProt ID:	O95544 , F5GXR5
Cytogenetics:	1p36.33
Synonyms:	dj283E3.1
Summary:	NADK catalyzes the transfer of a phosphate group from ATP to NAD to generate NADP, which in its reduced form acts as an electron donor for biosynthetic reactions (Lerner et al., 2001 [PubMed 11594753]).[supplied by OMIM, Mar 2008]
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism



[View online »](#)

Product images:

