

Product datasheet for TP720245XL

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

UROD (NM 000374) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human uroporphyrinogen decarboxylase (UROD)

Species: Human
Expression Host: E. coli

Expression cDNA Clone

Met1-Asn367

or AA Sequence:

Tag: N-His

Predicted MW: 43 kDa

Concentration: lot specific

Purity: >95% as determined by SDS-PAGE and Coomassie blue staining

Buffer: Provided lyophilized from a 0.2 μm filtered solution of 20 mM Tris-HCl, 150 mM NaCl

Endotoxin: < 0.1 EU per μg protein as determined by LAL test

Storage: Store at -80°C.

Stability: Stable for at least 3 months from date of receipt under proper storage and handling

conditions.

RefSeg: NP 000365

 Locus ID:
 7389

 UniProt ID:
 P06132

 Cytogenetics:
 1p34.1

 Synonyms:
 PCT; UPD

Summary: This gene encodes an enzyme in the heme biosynthetic pathway. This enzyme is responsible

for catalyzing the conversion of uroporphyrinogen to coproporphyrinogen through the removal of four carboxymethyl side chains. Mutations and deficiency in this enzyme are known to cause familial porphyria cutanea tarda and hepatoerythropoetic porphyria.

[provided by RefSeq, Aug 2010]

Protein Families: Druggable Genome

Protein Pathways: Porphyrin and chlorophyll metabolism





Product images:

