

Product datasheet for **TP720260L**

ASF1A (NM_014034) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ASF1 anti-silencing function 1 homolog A (<i>S. cerevisiae</i>) (ASF1A)
Species:	Human
Expression Host:	<i>E. coli</i>
Expression cDNA Clone or AA Sequence:	Met1-Met204
Tag:	N-T7&C-His
Predicted MW:	25.4 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_054753
Locus ID:	25842
UniProt ID:	Q9Y294
Cytogenetics:	6q22.31
Synonyms:	CGI-98; CIA; HSPC146



[View online »](#)

Summary:

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The protein is a key component of a histone donor complex that functions in nucleosome assembly. It interacts with histones H3 and H4, and functions together with a chromatin assembly factor during DNA replication and repair. [provided by RefSeq, Jul 2008]

Product images: