

Product datasheet for **TP720290XL**

CXADR (NM_001338) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coxsackie virus and adenovirus receptor (CXADR)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Leu20-Gly237
Tag:	C-His
Predicted MW:	25.1 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH ₂ O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	NP_001329
Locus ID:	1525
UniProt ID:	P78310
Cytogenetics:	21q21.1
Synonyms:	CAR; CAR4/6; HCAR
Summary:	The protein encoded by this gene is a type I membrane receptor for group B coxsackieviruses and subgroup C adenoviruses. Several transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene are found on chromosomes 15, 18, and 21. [provided by RefSeq, May 2011]



[View online »](#)

Protein Families: Druggable Genome, Stem cell - Pluripotency, Transmembrane

Protein Pathways: Viral myocarditis

Product images:

