

## Product datasheet for **TP720336**

### SLITL2 (VASN) (NM\_138440) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human vasorin (VASN)
Species:	Human
Expression Host:	HEK293
Expression cDNA Clone or AA Sequence:	Cys24-Pro575
Tag:	C-His
Predicted MW:	60.7 kDa
Concentration:	lot specific
Purity:	>95% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	Provided lyophilized from a 0.2 µm filtered solution of 20 mM Tris-HCl, 150 mM NaCl
Endotoxin:	< 0.1 EU per µg protein as determined by LAL test
Reconstitution Method:	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. Dissolve the lyophilized protein in ddH <sub>2</sub> O. It is not recommended to reconstitute a concentration less than 100 µg/ml. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
Storage:	Store at -80°C.
Stability:	Stable for at least 6 months from date of receipt under proper storage and handling conditions.
RefSeq:	<a href="#">NP_612449</a>
Locus ID:	114990
UniProt ID:	<a href="#">Q6EMK4</a>
Cytogenetics:	16p13.3
Synonyms:	SLITL2
Summary:	May act as an inhibitor of TGF-beta signaling.[UniProtKB/Swiss-Prot Function]
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)

**Product images:**

